

Tubing Lengths	Coiled Tubing		Stick Tubing
	Coated CP	Coated EP & Uncoated EP	Coated & Uncoated
1/8" X 0.020" SMLS 316L	900' (275M) Max	200' (60M) Max	20' (6M)
1/4" X 0.035" SMLS 316L	2,000' (610M) Max	600' (180M) Max	20' (6M)
3/8" X 0.035" SMLS 316L	1,300' (395M) Max	600' (180M) Max	20' (6M)



Common material and sizes shown. Call for other requirements.
 All tubing meets ASTM specifications: A269, A213-EAW, A1016, MR0175, EN 10204-3.1

SilcoTek. SilcoNert™ 2000 Coated Electropolished 316L Seamless Stainless Steel Tubing

- Ultimate Inertness
- Corrosive Environments
- Low Level Measurements

OD	Wall	Material	Stick or Coil	Part #	Description
1/8"	0.035	316L Seamless	Stick	020206LSEPS2	Electropolished 1/8" x 0.035 316L SMLS w/ SilcoNert 2000 in sticks
1/8"	0.035	316L Seamless	Coil	020206LSEPS2CR	Electropolished 1/8" x 0.035 316L SMLS w/ SilcoNert 2000 in coil
1/4"	0.035	316L Seamless	Stick	040356LSEPS2	Electropolished 1/4" x 0.035 316L SMLS w/ SilcoNert 2000 in sticks
1/4"	0.035	316L Seamless	Coil	040356LSEPS2CR	Electropolished 1/4" x 0.035 316L SMLS w/ SilcoNert 2000 in sticks

SilcoTek. SilcoNert™ 2000 Coated Chemically Polished & Passivated (CP) Tubing 316L Seamless Stainless Steel Tubing

- General Chromatography
- Inert Gas Delivery Systems
- Analytical Measurements

OD	Wall	Material	Stick or Coil	Part #	Description
1/8"	0.035	316L Seamless	Stick	020206LSBAS2	CP 1/8" x 0.035 316L SMLS w/ SilcoNert 2000 in sticks
1/8"	0.035	316L Seamless	Coil	020206LSBAS2CR	CP 1/8" x 0.035 316L SMLS w/ SilcoNert 2000 in coil
1/4"	0.035	316L Seamless	Stick	040356LSBAS2	CP 1/4" x 0.035 316L SMLS w/ SilcoNert 2000 in sticks
1/4"	0.035	316L Seamless	Coil	040356LSBAS2CR	CP 1/4" x 0.035 316L SMLS w/ SilcoNert 2000 in coil

Uncoated Electropolished 316L Seamless Stainless Steel Tubing

- Moisture Measurements
- Clean and High Purity
- Mirror Surface Finish
- Surface Roughness of 10 Ra µin

OD	Wall	Material	Stick or Coil	Part #	Description
1/8"	0.035	316L Seamless	Stick	020206LSEP	Electropolished 1/8" x 0.035 316L SMLS in sticks
1/8"	0.035	316L Seamless	Coil	020206LSEPCR	Electropolished 1/8" x 0.035 316L SMLS in coil
1/4"	0.035	316L Seamless	Stick	040356LSEP	Electropolished 1/4" x 0.035 316L SMLS in sticks
1/4"	0.035	316L Seamless	Coil	040356LSEPCR	Electropolished 1/4" x 0.035 316L SMLS in sticks

©2010 O'Brien Corporation • GLT-STOCKTUBE • 10 JUN 10



O'Brien Analytical & Cardinal Offices:

USA: 1910 Crystal Industrial Ct. • St. Louis, MO 63114 • Ph: +1 (314) 473-2121 • Fax: +1 (314) 473-2120

BELGIUM: Mallekotstraat 65 • B 2500 Lier Belgium • Ph: +32 (0)3 491 9875 • Fax: +32 (0)3 491 9876

CHINA: Suite 1604A, 16F Feidiao International Building • 1065 Zhao Jia Bang Road, Shanghai 200030 • Ph: +8621 51581635 • Fax: +8621 51581686

SINGAPORE: No. 14 Robinson Road • #13-00, Far East Finance Building • Singapore 048545 • Ph: +65 6387 8907 • Fax: +65 6387 9807

www.cardinaluhp.com • www.obrien-analytical.com